Phonak Insight●

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1. Windle reference in APA: Windle, R., Dillon, H., & Heinrich, A. (2023). A review of auditory processing and cognitive change during normal ageing, and the implications for setting hearing aids for older adults. *Frontiers in Neurology*, 14, Article 1122420. <https://doi.org/10.3389/fneur.2023.1122420>
2. United Nations Department of Economic and Social Affairs, Population Division (2022). *World Population Prospects 2022: Summary of Results*. UN DESA/POP/2022/TR/NO
3. World Health Organization. (n.d.). Ageing and health. Retrieved November 22, 2024, from <https://www.who.int/news-room/fact-sheets/detail/ageing-and-health>
4. Dryden, A., Allen, H. A., Henshaw, H., & Heinrich, A. (2017). The Association Between Cognitive Performance and Speech-in-Noise Perception for Adult Listeners: A Systematic Literature Review and Meta-Analysis. *Trends in Hearing*, 21, 2331216517744675. <https://doi.org/10.1177/2331216517744675>
5. Gatehouse, S., Naylor, G., & Elberling, C. (2006). Linear and nonlinear hearing aid fittings—2. Patterns of candidature. *International Journal of Audiology*, 45, 153-171. <https://doi.org/10.1080/14992020500429484>
6. Gatehouse, S., Naylor, G., & Elberling, C. (2006). Linear and nonlinear hearing aid fittings—1. Patterns of benefit. *International Journal of Audiology*, 45, 130-152. <https://doi.org/10.1080/14992020500429518>
7. Gatehouse, S., Naylor, G., & Elberling, C. (2003). Benefits from hearing aids in relation to the interaction between the user and the environment. *International Journal of Audiology*, 42(Supplement 1), S77-S85. <https://doi.org/10.3109/14992020309074627>
8. Derleth, P., Georganti, E., Latzel, M., Courtois, G., Hofbauer, M., Raether, J., & Kuehnel, V. (2021). Binaural Signal Processing in Hearing Aids. *Seminars in hearing*, 42(3), 206-223. <https://doi.org/10.1055/s-0041-1735176>
9. Woodward, J., Jansen, S., Kuehnel, V. (2020). Hearing inspired by nature: the new APD 20 fitting formula with adaptive compression by Phonak. *Phonak Insight*. Retrieved from www.phonak.com/evidence
10. Hassager, H. G., Wiinberg, A., and Dau, T. (2017a). "Effects of hearing-aid dynamic range compression on spatial perception in a reverberant environment," *J. Acoust. Soc. Am.* 141, 2556-2568
11. Kowalewski, B., Zaar, J., Fereczkowski, M., MacDonald, E., Strelcyk, O., May, T., & Dau, T. (2017). Effects of slow-and fast-acting compression on hearing impaired listeners' consonant-vowel identification in interrupted noise. In *International Symposium on Auditory and Audiological Research*. The Danavox Jubilee Foundation
12. Lopez-Poveda, E. A., Johannesen, P. T., Perez-González, P., Blanco, J. L., Kalluri, S., & Edwards, B. (2017). Predictors of hearing-aid outcomes. *Trends in hearing*, 21, 1-28
13. Ricketts, T. A. (2005). Directional hearing aids: then and now. *Journal of Rehabilitation Research and Development*, 42(4, Suppl 2), 133-144
14. Mueller, G. H., & Ricketts, T. A. (2000). Directional-microphone hearing aids: an update. *The Hearing Journal*, 53(5), 10-19
15. Bentler, R. A. (2005). Effectiveness of directional microphones and noise reduction schemes in hearing aids: A systematic review of the evidence. *Journal of the American Academy of Audiology*, 16(7), 473-484
16. Latzel, M., LeSimple, C., & Woodward, J. (2022). New implementation of directional beamforming configurations shows improved speech understanding and reduced listening effort *Phonak Field Study News*. Retrieved from www.phonak.com/evidence
17. Appleton, J. (2020). AutoSense OS 4.0 - significantly less listening effort and preferred for speech intelligibility. *Phonak Field Study News* retrieved from www.phonak.com/evidence
18. Latzel, M., Heeren, J and Lesimple, C. (2024). "Speech Enhancer reduces listening effort and fatigue." *Phonak Field Study News*. Retrieved from www.phonak.com/evidence
19. Raufer, S., Kohlhauer, P., Jehle, F., Kühnel, V., Preuss, M., Hobi, S. (2024). Spheric Speech Clarity proven to outperform three key competitors for clear speech in noise. *Phonak Field Study News*. Retrieved from www.phonak.com/evidence
20. Wright, A., Kuehnel, V., Keller, M., Seitz-Paquette, K., Latzel, M. (2024). "Spheric Speech Clarity applies DNN signal processing to significantly improve speech understanding from any direction and reduce the listening effort." *Phonak Field Study News* retrieved from <https://www.phonak.com/evidence>
21. Thibodeau L. M. (2020). Benefits in Speech Recognition in Noise with Remote Wireless Microphones in Group Settings. *Journal of the American Academy of Audiology*, 31(6), 404-411
22. Sonova proprietary research. (2019). Project ID#1299. Please contact marketinsight@phonak.com if you are interested in further information

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Stacey Rich,
Senior Manager Audiology Thought
Leadership.

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